

## Get more information from a single assay.

FusionPlex assays enable simultaneous detection of fusions, splice variants, key SNV/Indels, and relative RNA expression from the same input.

### Detection of multiple mutation types with FusionPlex



Cancer is a complex disease to study with multiple important biomarkers; therefore, a comprehensive discovery method is critical. The desire to optimize cost efficiency coupled with the need to preserve precious tissue samples is driving researchers to seek solutions that do more with less. An effective research assay must be multiplexed, user-friendly, and able to detect fusions, splice variants, and key SNV/Indels from the same, often limited, input in parallel. For fusion detection, RNA is superior to DNA due to the inherent limitations associated with amplification and mapping of intronic sequences.

Archer's FusionPlex assays offer a scalable and easy-to-use solution for RNA-based sequencing that meets the research needs of every lab. With Anchored Multiplex PCR (AMP™) chemistry at the core and an integrated bioinformatic platform, FusionPlex targeted NGS assays offer streamlined fusion detection along with SNV/Indel calling from sample to data. Choose from a suite of thoughtfully designed catalog panels or customize to meet your exact needs. FusionPlex assays can be used stand-alone or paired with a VariantPlex™ assay to generate a comprehensive genomic profile.

### FusionPlex NGS research assays offer advantages for customer workflows

	Simple, fast workflow	Known fusion detection	Novel fusion detection	Oncogenic isoform detection	Easy panel customization
FusionPlex® AMP™ chemistry (RNA)	✓	✓	✓	✓	✓
Hybrid capture (RNA)	✗ <sup>1</sup>	✓	✓	✓	✗
Hybrid capture (DNA)	✗ <sup>1</sup>	✓	! <sup>2</sup>	! <sup>4</sup>	✗
Opposing primer amplicon (RNA)	✓	✓	! <sup>3</sup>	! <sup>5</sup>	!


! Indicates that the methodology may not be the best choice for optimal results

1. Complex library preparation processes (Hsiao et al. 2019)
2. Introns can be challenging to target and map reads to resulting in coverage gaps (Benayed et al. 2019, Davies et al. 2018)
3. Breakpoint-spanning reads are not sequenced; expression imbalance-based analyses can result in uncertainty (Vendrell et al. 2017)
4. Isoform expression cannot be evaluated using DNA input
5. Relative mRNA abundance cannot be confirmed without molecular barcodes (MBCs)


## Engineered for reliability.

Patented Anchored Multiplex PCR (AMP) target enrichment chemistry provides robust detection of oncogenic drivers.

- Optimized for FFPE samples
- Known and novel fusion detection
- Molecular barcode (MBC)-driven error correction and unique molecule identification




## cDNA, DNA, or ctDNA fragments



## Surprisingly simple NGS workflow.

Designed with simplicity in mind, Archer makes it easy to implement NGS research assays in your lab.

- 1.5 day library prep with minimal hands-on time
- Input requirements as low as 10 ng RNA
- Lyophilized reagents to minimize error and eliminate the need for master mixes
- Single-use reactions reduce contamination
- Parallel workflow for RNA, DNA, and ctDNA analysis



## One powerful platform. Endless potential.

All FusionPlex assays are supported by an integrated analysis platform, Archer® Analysis, providing a turnkey solution to complex genomic problems.

- Simple and intuitive web-based interface
- In-line visualization for clear reporting
- Access to Quiver, our curated database of known fusion events
- Integrates LIMS data and third party providers
- Deploys securely to a cloud or local server



# Introducing Archer's newest FusionPlex panel

Now with expanded SNV/Indel content

FusionPlex Core Solid Tumor is optimized for expanded SNV/Indel detection

**FUSIONPlex**<sup>®</sup>  
Core Solid Tumor

AKT1	DDR2	FGFR2	IDH1	MAP3K8	NTRK3	PRKCA	TRIM11
ALK	DNAJB1	FGFR3	IDH2	MET	NUTM1	PRKCB	
AXL	EGFR	GNA11	KEAP1	MYB	PAX8	RAF1	
BRAF	ERBB2	GNAQ	KIT	MYBL1	PDGFRA	RET	
BRD3	ERBB4	GNAS	KRAS	NRAS	PIK3CA	ROS1	
BRD4	ERG	H3F3A	LTK	NRG1	POLD1	STK11	
CTNNB1	ESR1	HIST1H3B	MAP2K1	NTRK1	POLE	TMPRSS2	
CYSLTR2	FGFR1	HRAS	MAP3K3	NTRK2	PPARG	TP53	


SNV/Indel targets

## Sensitive fusion detection from samples with low tumor content

Fusion and oncogenic isoform detection in FFPE reference material at low input masses

- $\geq 95\%$  of expected fusions & oncogenic isoforms detected from 50 ng of 10% fusion-positive input.

Reference materials used: (1) *Seraseq*<sup>®</sup> FFPE Fusion RNA Reference Material v4 (0710-0496), (Background) *Seraseq*<sup>®</sup> FFPE WT (DNA/RNA) Reference Material (0710-0137);  
n = 3 per input mass, per dilution.



## Sensitive SNV/Indel detection from RNA with high concordance to DNA

Libraries were generated from FFPE-derived TNA using RNA-based FusionPlex and DNA-based VariantPlex assays.

93 SNV/Indels with  $\geq 5\%$  VAF were identified by the DNA-based assay.

- RNA:DNA concordance was 85.7% for 32 poor quality samples and 100% for 47 moderate quality samples.


Gene	42 variants detected in 32/42 poor quality* samples								51 variants detected in 47/51 moderate quality* samples							
	TOTAL	EGFR	ERBB2	FGFR3*	KRAS	MET	NRAS	PIK3CA	TOTAL	BRAF	EGFR	ESR1	MET	PDGFRA	RAF1	RET
# DNA variant (VariantPlex)	42	6	1	5	21	5	1	3	51	4	6	1	4	34	1	1
# RNA variant (FusionPlex)	36	6	1	1*	21	4	1	2	51	4	6	1	4	34	1	1
Concordance (%)	85.7%	100%	100%	20%	100%	80%	100%	67%	100%	100%	100%	100%	100%	100%	100%	100%

\*Quality assessments for RNA are based on in-line PreSeq<sup>®</sup> RNA QC scores. Overall quality assessments are based on library complexity (total unique fragments).

## Fusion and SNV/Indel detection from the same TNA input


FusionPlex research assays conserve precious tissue by providing a comprehensive genomic profile from a single TNA input. SNVs that confer resistance to tyrosine kinase inhibitors were identified in *NTRK* fusion positive FFPE tissue.

<i>NTRK</i> fusion	<i>NTRK</i> resistance mutation (SNV)	Expressed SNV AF (%)
<i>ETV6:NTRK3</i>	<i>NTRK3</i> G623R (chr15:88476265 C>T)	0.39%
<i>ETV6:NTRK3</i>	<i>NTRK3</i> G623R (chr15:88476265 C>T)	1.75%
<i>ETV6:NTRK3</i>	<i>NTRK3</i> G623R (chr15:88476265 C>T)	0.39%



## Thoughtfully designed catalog assays plus easy customization to meet your research needs

Compatible with Illumina® and Ion Torrent™ sequencers

Assays	Illumina catalog #	Ion Torrent catalog #	<b>Custom:</b> Design your own panel or add to an existing product. It's your panel, your way. 
FusionPlex Core Solid Tumor	AB0145	AB0146	
FusionPlex Lung v2	AB0135	AB0136	
FusionPlex Pan Solid Tumor v2	AB0137	AB0138	
FusionPlex Sarcoma v2	AB0133	AB0134	
FusionPlex Myeloid	AB0015	AB0016	
FusionPlex Pan Heme	AB0017	AB0018	
FusionPlex Heme v2	AB0011	AB0012	
FusionPlex ALL	AB0013	AB0014	
FusionPlex Lymphoma	AB0019	AB0020	

Learn more at [www.archerdx.com/solid-tumor-research](http://www.archerdx.com/solid-tumor-research) or email us at [sales@archerdx.com](mailto:sales@archerdx.com)

**For research use only. Not for use in diagnostic procedures.** Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

© 2023 Integrated DNA Technologies, Inc. All rights reserved. Trademarks contained herein are the property of Integrated DNA Technologies, Inc. or their respective owners, and may be registered in the USA and/or other jurisdictions. For specific trademark and licensing information, see [www.idtdna.com/trademarks](http://www.idtdna.com/trademarks). RUO23-1815\_001.2 04/14